

Pathology Division

Section: Anatomic Pathology / Technical Procedures		Number: PRC-PAT(IHC)-500
Title: Transportation of Direct Immunofluorescence Slides		Page 1 of 1
Issuing Authority (s): B.Dyer, Pathology Manager Signature: _____ Date: Sept. 25, 2008		Author: B.Rowe
Issue Date: Sept. 26, 2008	Date Effective: Sept. 26, 2008	Revision Dates:

Purpose

This procedure provides instructions for transportation of Direct Immunofluorescence slides from HSC Pathology lab to St. Clare’s Pathology lab.

Materials

Supplies

Equipment

- Blue transportation bag

Procedure

Follow the activities below to send Direct Immunofluorescence slides to the Pathologist.

1. Photocopy Immunofluorescence Request form. Place copy in ‘Renal Bx’ binder. Place original form in the folder with the Immunofluorescence slides.
2. When the H&E and special stains are completed the slides MUST be delivered to the IHC Department.
3. Remove Immunofluorescence slides from IHC fridge.
4. Place all slides, requisition and worksheets in a blue transportation bag.
5. Call cab and bring blue transportation bag to Lab office at HSC.
6. Record information on Shipping bag/Sending log in Lab office.
7. Call St. Clare’s Pathology Lab to inform them that the slides are being transported.
8. Document the time and person called under ‘Spec Comments’ in Meditech.

Documents/Forms/Records

Immunofluorescence Request form