

P.A.P.

Day I

Prepare Slides:

- Cut sections @ 3 microns on clear glass slide.
- Place in 58-60 degree C. incubator overnight

Day II

- Place 200 ml 1/10 PBS in a 37 degree C. incubator.
- Weigh out 0.2 gm trypsin.
- Make up albuminized 1/10 PBS:

- 1 ml bovine albumin per 100 ml 1/10 PBS (This solution is used to prepare all antibody dilutions and for washing slides).

Procedure

(1) Transfer slides directly from 58-60 degree C. incubator to:

1. xylene 10 min.
2. xylene 10 min.
3. xylene 10 min.
4. xylene 10 min.
5. 95% alcohol 5 min.
6. Cold Tap H₂O 5-10 min.
7. 37 degree Tap H₂O 3-5 min.

(2) Add 0.2 gm trypsin to 200 ml 1/10 PBS in 37 degree C. incubator.

(3) Transfer slides to trypsin solution for 6 minutes at 37 degrees C. (Discard after use).

(4) Rinse in cold tap H₂O for 5 min.

Inhibition of Endogenous Peroxidase

(5) Transfer slides to freshly prepared 200 ml of methanol + 2 ml 30% H₂O₂ for 30 minutes

(6) Rinse in cold tap water for 5 minutes

(7) 95% alcohol 5 min.

- 2 -

- (8) Abs. alcohol 5 min.
- (9) Chloroform 5 min.
- (10) Acetone 5 min.
- (11) Cold Tap H₂O 5 min.
- (12) Albuminized PBS 5 min. minimum

Primary Antibody

- (13) Prepare dilutions of primary antibody using albuminized PBS 140 - 200 ul per slide.
- (14) Prepare slides one by one by drying outside and leaving a film of the albuminized PBS around the section (NO AIR BUBBLES). (DO NOT LET THE SECTIONS DRY!).
- (15) Cover section with diluted primary antibody (one slide at a time) for 30 minutes in a humidification chamber.
- (16) Wash each slide carefully with 1/10 PBS (albuminized) using a squeeze bottle.
- (17) Transfer to a staining dish of the albuminized PBS.

*Secondary Antibody

- (18) Prepare dilution as follows:
 for 20 slides
 PBS 1/10 - 2700 ul
 Normal Human - 300 ul
 Serum

Shake well and discard 200 ul of this solution. Add 200 ul peroxidase conjugated rabbit immunoglobulines to mouse immunoglobulins. Mix well.

- (19) Prepare slides as for primary antibodies and add secondary antibodies for 30 minutes.
- (20) Wash with albuminized PBS and return to staining dish of PBS.

Zymed

- 3 -

Teritary Antibody

- (21) Prepare dilutions as for secondary antibody using peroxidase conjugated swine immunoglobulins to rabbit immunoglobulins.
- (22) Prepare slides as before and add antibody for 30 minutes.
- * (23) Wash well with PBS and return to staining dish of PBS.

Revelation of Peroxidase Staining:

- (24) Prepare DAB:

Dissolve 6 mg of DAB (diaminobenzidine) for 6 slides in 10 ml Tris-HCL (0.05M). Add 3 ul of 30% H₂O₂ just before use. Keep in dark. (If larger quantities are needed, use several tubes ready for H₂O₂.)

H
10 ml

H
10 ml

H
10 ml

- (25) Arrange slides in humidification chamber. Put DAB on each slide for 3 to 5 minutes. (Do not attempt to reveal too many slides at one time - maximum 15 slides).

- (26) Rinse with old tap H₂O to stop reaction.

- (27) Return to a staining jar of water.

Counterstain

- (28) Haematoxylin for 10 to 20 seconds

- (29) Rinse in cold tap H₂O.

- (30) 95% alcohol

- (31) Absolute alcohol

- (32) Acetone

- (33) Xylene

Mount with Permount